



db

ф

Prepared to be done in the presence of the shower accessory.

B - Tube outlet for chromotherapy feeding (Ø20mm).
Prepare to be made in the presence of the 30x30 shower head with chromotherapy;

 $_{O}^{*}$ - Attacks A and B not to be considered in the presence of the G-Steam column.

ATTENTION:

Gemini Box with RIGENERA column or G-STEAM see its technical datasheet for attack predisposition.

- * The floor units for the connection of hot and cold water vary with the height of the shower tray.

WARNING: In case of non-observance, even partial, of safety conditions, both the guarantee and the liability of the manufacturer for any direct or indirect damage caused by the product. The loss of warranty, in the aforementioned hypothesis, will have the effect of replacing only pieces that will be recognized defective, without covering labor costs resulting from the removal of the product.

GRUPPO GEROMIN srl reserves the right to make changes to the product at any time without giving any notice